



The performance of FVIII and FIX measurement in ECAT surveys

Piet Meijer
ECAT Foundation
Leiden
The Netherlands



Introduction

The measurement of Factor VIII (FVIII) and Factor IX (FIX) clotting activity is used for both the diagnosis of haemophilia A or B and the monitoring of treatment and therefore requires precise laboratory measurement. We evaluated the performance of FVIII and FIX measurement in surveys of the ECAT external quality assessment programme.




Methods

- **Test results of FVIII and FIX of 175 - 215 participants were evaluated for the period 2010 – 2012.**
- **The between-laboratory variation (BCV) was assessed at different Factor VIII levels.**
- **The long-term analytical coefficient of variation (LCVa) was assessed using a linear regression model (P. Meijer *et al.* Clin Chem 2002;48:1011-15). The LCVa is a measure for the long-term analytical performance of a laboratory.**




Results FVIII + FIX






Results FVIII






Results FVIII





Results FVIII





Results FVIII

Severe Haemophilia A

N = 201


98x < [value] : 16 x < [0.1 – 0.9]
69 x < 1
13 x < [2 – 12] (13%)

103 numerical results

Robust mean : 0.95%
Median : 1.00%
SD : 0.75%
Range : 0 – 41%

Result > 1% : 22 (21%)

Overall: 17% reported a result > 1% FVIII





Discussion

Table 1. FVIII:C activity assay with one-stage method of five samples from severe haemophilia A patients (1, 5, 7, 8 and 10) and five samples from moderate haemophilia A patients (2, 3, 4, 6 and 9). Eleven different aPTT reagents were used in the FVIII:C assays.


Sample	Mutation gene type	Reagents										
		1	2	3	4	5	6	7	8	9	10	11
1	inv. Intron 1	1.4	0.4	0.4	0.3	0.9	0.1	0.3	0.8	1.4	0.3	0.5
2	ex10 Tyr 476 stop C → G	6.4	4.1	5.1	4.3	5.4	3.8	4.6	4.6	6.8	5.2	5.8
3	ex7 Hs256Arg A → G	2.5	1.4	1.5	1.4	2.1	1.2	1.3	1.9	3.6	1.7	1.8
4	ex14 Asp121461u C → G	5.4	3.9	4.5	3.7	5.0	3.8	4.1	4.6	8.3	5.0	4.9
5	inv. Intron 22	1.3	0.4	0.4	0.4	0.9	0.1	0.2	0.8	2.0	0.4	0.4
6	ex23 Arg2163Cys C → T	1.8	0.8	0.9	0.9	1.3	0.5	0.8	1.5	3.0	1.0	1.0
7	ex14 Arg795stop C → T	nd	0.4	0.7	0.3	0.6	0.1	0.3	0.9	1.3	nd	0.5
8	inv. Intron 1	1.7	1.0	1.0	0.8	0.9	0.5	0.6	0.9	1.8	0.7	0.9
9	ex7 His256Arg A → G	2.8	1.7	2.0	1.3	3.0	1.6	1.9	2.9	1.9	2.1	2.6
10	inv. Intron 22+/-	2.1	1.1	0.4	1.0	1.9	1.0	1.3	2.0	1.6	1.2	1.6

The results of the assays are expressed as IU dL⁻¹.

Reagents: 1, STA Cephascreen (Stago); 2, STA APTT Kaolin (Stago); 3, PTT-LA (Stago); 4, PTT a (Stago); 5, PTT Reagent (Stago); 6, Synthasil (HemosIL); 7, APTT-SP (HemosIL); 8, Dade[®] Actin[®] FSL Activated PTT Reagent (Dade Behring); 9, Dade[®] Actin[®] FS Activated PTT Reagent (Dade Behring); 10, Pathromtin SL (Dade Behring); 11, Platelin[®] LS (Biomerieux); NO, Not Determined.




Results FIX





Results FIX





Results FIX

Severe Haemophilia B

<u>Factor IX (plasma 10.60)</u>	n	mean	CV	range
total group	154	1.1	61.1%	0.1 – 4.3
One-stage Clotting Assay	152	1.1	61.1%	0.1 – 4.3


<u>Activator</u>	n	mean	CV	Range
HemosIL APTT-SP liq. sil. (IL)	10	1.0	18.2%	0.5 – 1.3
HemosIL Synthasil (IL)	20	1.0	47.6%	0.3 – 2.0
Actin FS (Siemens)	16	1.1	45.9%	0.8 – 3.0
Actin FSL (Siemens)	16	1.0	42.3%	0.2 – 2.0
Pathromtin SL (Siemens)	20	1.1	103%	0.1 – 4.3
PTT Automate./STA APTT (Stago)	24	1.1	56.9%	0.3 – 3.0
Ceph./Kaolin - CKPrest (Stago)	15	1.0	58.3%	0.1 – 3.0

18x < [value]

Numerical results : 29 x result > 1% (19%)



Long-term evaluation



X = consensus value ; \bar{X} = mean value for X .

s_x = standard error of X

Y = laboratory value ; \bar{Y} = mean value for Y .

b = slope

$s_{y|x}$ = variability of the regression line, which is calculated based on the least-square method.

n = number of laboratory results

IMPRECISION

$$LCV_a = \frac{s_{y|x}/b}{\bar{X}} \cdot 100\%$$

BIAS


$$B = \frac{\sqrt{\frac{n-1}{n} \cdot b-1^2 \cdot s_x^2 + \bar{Y} - \bar{X}^2}}{\bar{X}} \cdot 100\%$$

TOTAL ERROR

$$TE = \frac{\sqrt{s_{y|x}^2 + b-1^2 \cdot s_x^2 + \bar{Y} - \bar{X}^2}}{\bar{X}} \cdot 100\%$$



Long-term evaluation FVIII



Descriptive statistics

N : 187
Mean : 15.7%
Median : 13.3%
Range : 5.1 – 39.9%
(2.5 – 97.5 CI)

Performance goals
based on the biological
variation

Diagnosis Monitoring

Factor
VIII

17.5

8.3

% Pass


71%

16%



Long-term evaluation FVIII


Chromogenic



Long-term CVanalytical 5.8%

Bias 12.1%

Clotting



Long-term CVanalytical 5.6%

Bias 0.5%



SIX-SIGMA

$$\text{Sigma} = (\text{TEa} - \text{Bias}) / \text{CV}$$

$$\text{TEa} < 1.65 \times (0.5 \times \text{CV}_{\text{ws}}) + 0,25 \times (\text{CV}_{\text{ws}}^2 + \text{CV}_{\text{bs}}^2)^{1/2}$$

Factor VIII

TEa = 13.3% (Westgard database – minimum specs)

Chromogenic

Bias = 12.1%

CV = 5.8%

Sigma = 0.21

Clotting

Bias = 0.5%

CV = 5.6%

Sigma = 2.29



Conclusions

- **The between-laboratory variation for FVIII and FIX is comparable.**
- **Differences between reagents are small.**
- **The accuracy of measurement in the low range needs improvement.**
- **The long-term performance evaluation may help laboratories in their quality management.**
- **Especially for monitoring treatment a more stable performance over time is needed.**